



**Power to Control Bleeding
Power to Heal**



Sealk™

**Advanced
Hemostatic
Dressing**



For Ear, Nose and Throat

100% Natural Plant Based Water Soluble

Etherified Carboxymethyl Cellulose

**Power to Control Bleeding.
Power to HEAL.**

For Otolaryngology

For use to promote healing in these procedures:

- Control of Epistaxis
- Tonsillectomy and/or Adenoidectomy
- Septoplasty
- Functional Endoscopic Sinus Surgery (FESS)
- Functional Rhinoplasty
- Thyroid¹

¹BloodSTOP iX has European Union CE Certification for implant surgical and topical use, including when used with anticoagulant medications. BloodSTOP iX currently has USFDA clearance for topical wound hemostasis.



Sealk is an all-natural, plant-based, biocompatible, water soluble etherified carboxymethyl cellulose matrix that utilizes patented technology to achieve superior hemostasis. It contains no animal, human, or added chemical components.

64% of post-tonsillectomy patients reported less pain on the side treated with **Sealk**¹

Average reduction in pain is¹ **>40%**

Pain relief seems to be consistent¹

DAY 1-10

¹Bunting, E. Hawley, B. Therapeutic Effects of Sealk on Post Tonsillectomy Wound Sites, 2023, Clinical Trial, 50 Pts.

Contact us for more information:

- 📍 1356 N Emporia Ave
Wichita, KS 67214
- 📞 316-251-0920
- 📠 fax: 316-251-0950
- ✉ purchase@thomasjanemedical.com

Sealkproducts.com

Distributor: **Thomas Jane MEDICAL**



What is Sealk?

Sealk is a clinically proven, research-based hemostatic dressing and wound healing matrix that provides rapid hemorrhage control, dramatically reducing blood loss. It is:

- All-natural, plant-based
- Biocompatible
- Water soluble
- Activates growth factors, promoting tissue regeneration and angiogenesis²
- Decreases inflammatory cytokines, reducing the infection and inflammation risk²
- Promotes fast wound healing³

²Peng D, B. Reed-Maldonado A, Banie L, Wang G, Lin G, F. Lue T. Carboxymethylcellulose Activates Dermal Cells and Adipose-Derived Stem Cells Through Wnt/ β -catenin Pathway. *J Surg Res (Houst)*. 2021;04(01). doi:10.26502/jsr.10020117
³Ju S, Wang K, Qiao L et al. Application of BloodSTOP iX Wound Heal Nanocellulose Matrix for Burn Wound Care. *J Surg Res (Houst)*. 2021;04(01). doi:10.26502/jsr.10020105



Sealk adhered to RIGHT tonsil bed post-tonsillectomy.



Safe, Natural & Effective
for use on a Variety of Patients

It's Easy to Use!

1 Keep Dry



2 Cut to Size as Needed



3 Place on Wound



4 Irrigate as needed to seal



**VIEW THE DEMO
& LEARN MORE**



Features & Benefits

- Rapid hemostasis
- Translucent gel for easy wound monitoring
- The gel sticks, seals & protects wounds
- Prevents contamination, conforms to and adheres to wound sites
- pH neutral: ~7.2, biocompatible
- Easily applied and removed without irritation
- Provides a moist wound environment to promote healing
- Leaves no residue behind
- Highly cost effective
- Reduces pain upon application

Give your Patients Peace of Mind

“

TESTIMONIAL

I am no longer prescribing narcotics for pediatric T&A after implementing Sealk. It works so well that patients and parents are requesting it now!

Dr. K. Smothers, ENT

”